



TechniCare
Improving Healthcare



TOBA 2 Shower chair

Newest generation shower chair



“Designed for common use to high-complex care”

This is the TOBA 2. By using the newest systems and technics, the TOBA 2 can be used for both common use and client-specific solutions. By having a large number of accessories, the TOBA 2 can be adjusted to wishes and needs of the client.

The TOBA 2 is the newest generation shower-toilet chair with an electrical

in height and tilt adjustment to secure safe and ergonomic washing and caring by the professional. During the height adjustment the shower chair is automatically tilting to secure the best seating position.

We have developed the TOBA 2 for common use and thanks to diverse

accessories the TOBA 2 can be adjusted to clients' needs. On the back and seating side of the TOBA 2 diverse paddings can be installed.

With foldable foot rest and arm rests with front closure we deliver safety and easy access. Because of this, the TOBA 2 has access from both the side and front.



Advantages

Electrical height adjustment with independent tilt function

Central brake system

Foldable foot plate with transparent antislip coating

Adjustable foot plate to the floor

In height adjustable back section for maximum comfort.

Removable and foldable arm rests

Tilt section up to 15°

Easy remote control

Newest electrical components

Improving Healthcare

We believe in a better healthcare. TechniCare has the mission to improve the healthcare sector by technical innovation. We deliver high quality products, developed in association with the healthcare sector. We are a Dutch company with our own developing and producing facilities. Please find more information at www.technicare.nl



Nikkelstraat 24
2984 AM Ridderkerk
The Netherlands

T +31 (0)85 087 04 21
info@technicare.nl

www.technicare.nl