



**TechniCare**  
Improving Healthcare



# Standing aid Jill

Transfer and rehabilitation



## Standing aid Jill

The standing aid is designed for clients in care facilities who can stand up independently but have difficulty walking or are unable to walk (mobility classes A, B, and C). The standing aid enables caregivers to assist clients from a seated position in a safe and ergonomic way, helping them to stand and be moved over short distances. This includes transfers, for example, from the bed to a toilet or (wheel)chair.



### ✓ Horizontal foot spread

So that you can use the aid with any wheelchair, shower chair, or care chair.

### ✓ Rehabilitation function

Because the patient can pull themselves up as much as possible using the pull-up supports, the aid is directly used as a rehabilitation tool.

### ✓ Sit-to-stand function

The seat cushions are foldable and collapsible. They are positioned in such a way that the patient can stand for an extended period, enhancing the functionality.

## Electric leg spreading

The width of the electrically adjustable leg spreading can be easily adjusted using the hand control to accommodate different situations, allowing the JILL to be positioned close to the client at all times. This promotes a smooth transfer.

Once the feet are placed on the sturdy anti-slip footplate, the client can pull themselves up using the comfortable handles located at the top of the standing aid and stand up.

The JILL actively encourages the client to participate in the transfer using their own strength, thereby increasing their independence.

During use, the two swivel wheels can be locked.



### Specifications of the standing aid:

Minimum total width:	700 MM
Maximum total width:	1000 MM
Seat height:	857 MM
Seat width:	535 MM
Total length:	925 MM
Total height:	1208 MM
Inner width of leg spread:	Min. 585 MM, Max. 885 MM

# Improving Healthcare

TechniCare aims to improve healthcare by providing technical innovation. TechniCare delivers high-quality products developed in collaboration with the end user. We are a Dutch company with our own development and manufacturing facilities. For more information:

[www.technicare.nl](http://www.technicare.nl)



**TechniCare**  
Improving Healthcare

Nikkelstraat 24  
2984 AM Ridderkerk

T +31 (0)180 499 276  
info@technicare.nl

[www.technicare.nl](http://www.technicare.nl)